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## Literatúra

- Ayala-Cabrera J.F., Contreras L., Moyano E., Santos F.J. (2020) A novel methodology for the determination of neutral perfluoroalkyl and polyfluoroalkyl substances in water by gas chromatography-atmospheric pressure photoionisation-high resolution mass spectrometry. *Anal. Chim. Acta* DOI: 10.1016/j.aca.2019.12.004.
- Dauchy, X. Bioteux V., Back C., Colin A., Hemard J., Rosin C., Munox J., (2017) Mass flows and fate of per- and polyfluoroalkyl substances (PFASs) in the wastewater treatment plant of a fluorochemical manufacturing facility *Sci. Total Environ.* 576 549-558.
- Dimzon I.K., Wsterveld J., Gremmel C., Fromel T., Knepper T.P., de Voogt P. (2017) Sampling and simultaneous determination of volatile per- and polyfluoroalkyl substances in wastewater treatment plant air and water *Anal Bioanal Chem* 409: 1395-1404.
- Favreau, P.; Poncioni-Rothlisberger, C.; Place, B. J.; Bouchex- Bellomie, H.; Weber, A.; Tremp, J.; Field, J. A.; Kohler, M. Multianalyte Profiling of Per- and Polyfluoroalkyl Substances (PFASs) in Liquid Commercial Products. *Chemosphere* 2017, 171, 491–501.

- Higgins, C.; Field, J.; Deeb, R.; Conder, J. FAQs Regarding PFASs Associated with AFFF Use at U.S. Military Sites; Environmental Security Technology Certification Program Alexandria United States, 2017.

- Herzke, D.; Olsson, E.; Posner, S. Perfluoroalkyl and Polyfluoroalkyl Substances (PFASs) in Consumer Products in Norway – A Pilot Study. *Chemosphere* 2012, 88, 980–987.

- Kim, M. H.; Wang, N.; McDonald, T.; Chu, K.-H. Biodefluorination and Biotransformation of Fluorotelomer Alcohols by Two Alkane-Degrading Pseudomonas Strains. *Biotechnol. Bioeng.* 2012, 109, 3041–3048.

- Ma H., Peng H., Chen H., Shang W., Zheng X., Yang M., Zhang Y., (2022) Long-term trends of fluorotelomer alcohols in a wastewater treatment plant impacted by textile manufacturing industry, *Chemosphere*, Volume 299.

- Portolés T., Rosales L.E., Sancho J.V., Santos J., Moyano E., (2015) Gas chromatography-tandem mass spectrometry with atmospheric pressure chemical ionization for fluorotelomer alcohols and perfluorinated sulfonamides determination, *Journal of Chromatography A*, Volume 1413, 2015, 107-116.

- Titley I.A., Florentino B., Cruz D., Barlaz M., Field J.A. (2023) Neutral Per- and Polyfluoroalkyl Substances in In-situ Landfill Gas by Thermal Desorption-Gas Chromatography-Mass Spectrometry *Environ. Sci. Technol. Lett.* 2023, 10, 3, 214-22.

- Wang, N.; Szostek, B.; Buck, R. C.; Folsom, P. W.; Sulecki, L. M.; Capka, V.; Berti, W. R.; Gannon, J. T. Fluorotelomer Alcohol Biodegradation Direct Evidence That Perfluorinated Carbon Chains Breakdown. *Environ. Sci. Technol.* 2005, 39, 7516–7528.

- Yan P.F., Dong S, Manz K.E., Liu C., Woodcock M.J., Mezzari M.P., Abriola L.M., Pennell K.D., Cápiro N.L. Biotransformation of 8:2 Fluorotelomer Alcohol in Soil from Aqueous Film-Forming Foams (AFFFs)-Impacted Sites under Nitrate-, Sulfate-, and Iron-Reducing Conditions. *Environ Sci Technol.* 2022 Oct 4;56(19):13728-13739. doi: 10.1021/acs.est.2c03669.

Kontaktujte našich expertov